

# Solvent Recovery Unit

Model: GAS510C

- First edition -

● Thank you for purchasing "Solvent Recovery Unit GAS510C" of Yamato Scientific Co.,Ltd.

● To use this unit properly, read this "Instruction Manual" thoroughly before using this unit.  
Keep this instruction manual around this unit for referring at anytime.

 : WARNING!

Carefully read and thoroughly understand the important warning items described in this manual before using this unit.

**Yamato Scientific Co.,Ltd.**

<b>1. Safety precautions</b> .....	<b>1</b>
Explanation of pictograms .....	1
List of symbols.....	2
Warning · Cautions.....	3
<b>2. Before using this unit</b> .....	<b>5</b>
Precautions when installing the unit .....	5
<b>3. Names and functions of parts</b> .....	<b>9</b>
Main unit.....	9
GAS510C Internal mechanism .....	10
Operation interface overview.....	11
Online operation interface overview (ADL312SC as an example) .....	12
Description of switch and indicator lamp in the interface.....	13
Value display and input description in the interface.....	14
<b>4. Operating procedures</b> .....	<b>15</b>
Preparations .....	15
Preparations before operation (ADL312SC+GAS510C) .....	16
Preparation before operation (GB211C+GAS510C) .....	18
Preparation before operation (DL411C+GAS510C).....	20
Operating method.....	22
Gas source selection .....	29
The relation between blower and temperature/dry air amount (reference) .....	30
Cooling recovery capacity (reference).....	31
<b>5. Handling Precautions</b> .....	<b>35</b>
About applicable organic solvents .....	36
Corrosion resistance table .....	38
Precautions during operation.....	39
<b>6. Maintenance Method</b> .....	<b>41</b>
Daily Inspection and Maintenance.....	41
<b>7. Long storage and disposal</b> .....	<b>44</b>
When not using this unit for long term / When disposing .....	44
Matters to consider when disposing of the unit .....	44
<b>8. When a trouble occurs</b> .....	<b>45</b>
Safety unit and error indications .....	45
Trouble Shooting .....	48
<b>9. After Service and Warranty</b> .....	<b>50</b>
When requesting a repair .....	50
<b>10. Specifications</b> .....	<b>51</b>
<b>11. Wiring diagram</b> .....	<b>52</b>
Piping system diagram .....	53
List of materials .....	55
<b>13. Action principle</b> .....	<b>56</b>
<b>14. Replacement parts table</b> .....	<b>57</b>
<b>15. List of Dangerous Substances</b> .....	<b>59</b>

# 1. Safety precautions

## Explanation of pictograms


### About pictograms


A variety of pictograms are indicated in this operating instruction and on products for safe operation. Possible results from improper operation ignoring them are as follows.

Be sure to fully understand the descriptions below before proceeding to the text.

---

---

 **Warning** Indicates a situation which may result in death or serious injury (Note 1.)

 **Caution** Indicates a situation which may result in minor injury (Note 2) and property damages (Note 3.)

---


---


(Note 1) Serious injury means a wound, an electrical shock, a bone fracture or intoxication that may leave after effects or require hospitalization or outpatient visits for a long time.


(Note 2) Minor injury means a wound or an electrical shock that does not require hospitalization or outpatient visits for a long time.

(Note 3) Property damage means damage to facilities, devices and buildings or other properties.

### Meanings of pictograms

 This pictogram indicates a matter that encourages the user to adhere to warning ("caution" included).  
Specific description of warning is indicated near this pictogram.

 This pictogram indicates prohibitions  
Specific prohibition is indicated near this pictogram.

 This pictogram indicates matters that the user must perform  
Specific instruction is indicated near this pictogram.

# 1. Safety precautions

## List of symbols

### Warning



General warnings



Danger!: High voltage



Danger!: High temperature



Danger!: Moving part



Danger!: Hazard of explosion

### Caution



General cautions



Electrical shock!



Burning!



Caution for no liquid heating!



Caution for water leak!



For water only



Poisonous material

### Prohibitions



General bans



Fire ban



Do not disassemble



Do not touch

### Compulsions




General compulsions



Connect ground wire



Install levelly



Pull out the power plug



Regular inspection

# 1. Safety precautions

## Warning - Cautions

### Warning



#### **Do not use this unit in an area where there is flammable or explosive gas**

Never use this unit in an area where there is flammable or explosive gas. This unit is not explosion-proof. Always try to assure sufficient ventilation in the room and take extreme care so that the atmosphere will not reach the explosive limit concentration. See "15. List of Dangerous Substances" on P. 59 for explosive or flammable gases.



#### **Always ground this unit**

Always ground this unit on the power equipment side in order to avoid electrical shock due to a power surge.



#### **Apply the source of rated power or more**

Be sure to apply the source of rated power or more. Applying non-rated voltage or non-rated power supply may cause the fire or electric shock.



#### **Prohibition of use for error**

If a smoke or abnormal smell may be occurred, turn off the power switch of the main unit immediately, and turn off the original power source, and finally contact to either the dealer you purchased this unit or our sales office. Leaving the failure may cause the fire or electric shock. Since the repairing of this unit is dangerous for non-specified service person, never repair the unit by the customer himself.



#### **Do not use the power cord if it is bundled or tangled**

Do not use the power cord if it is bundled or tangled. If it is used in this manner, it can overheat and fire may be caused.



#### **Do not damage power cord**

Do not damage power cord by bending, pulling, or twisting forcedly. It may cause the fire or electric shock. Besides, operating the unit with the something put on the cord may cause overheat, and result in fire.



#### **Do not use flammable or explosive substances.**

Do not use flammable or explosive substances or substances containing flammable and explosive materials, which may cause explosion and fire.  
Please refer to P.59 [15. List of Dangerous Substances].



#### **Never try to touch a hot part.**

Some parts of the unit are hot during and immediately after operation. Take special care for possible burning.



#### **Never try to disassemble or alter the unit.**

Never try to disassemble or alter the unit. A malfunction, a fire or an electrical shock may result.

# 1. Safety precautions

Warning - Cautions

## Caution

### **During a thunder storm**

During a thunderstorm, turn off the power key immediately, then turn off the circuit breaker and the main power. If this procedure is not followed, fire or electrical shock may be caused.

### **After power outage**

When power is shut off due to turning the ELB "OFF" or a power failure during operation (while the blower is in operation or liquid is being sent), after the power is restored, the system will return to the standby state. Please firstly observe the inlet temperature of the spray dryer and GAS oxygen concentration. If the inlet temperature is higher than 100°C and the GAS oxygen concentration is higher than 8%, do not click the N<sub>2</sub> introducing device, which may cause explosion. When the inlet temperature drops below 100°C, it's able to click the N<sub>2</sub> introducing device.

### **Do not perform unattended operation during activating the unit**

Do not perform unattended operation during activating the unit. Since the unit is in idling status and the nozzle is blocked or after the operation using sample, the temperature around outlet is increased and the remaining sample is flown from the sample tube disconnected from the unit, and this failure may cause the indeterminism accident.

### **Any people other than the qualified personnel shall never attempt to operate the unit.**

Take sufficient care for the control of the unit so that any people other than the qualified personnel shall never have a chance to operate the unit.

### **Always monitor and check the oxygen concentration.**

Always monitor the oxygen concentration in the unit to assure safety.

### **Take care when opening the unit.**

Be sure to confirm that the oxygen concentration has returned to 21% and avoid putting your face close to the exhaust port carelessly when opening the unit.

### **Notes when using solvents.**

#### **The unit has been designed to use ethanol.**

Service lives of a filter element or packing may be influenced depending on the type of solvent used. When gas leakage or other trouble occurs inside the unit, immediately replace the defective part with a new one. Check whether a solvent can be used or not, please refer to P. 36 "About applicable organic solvents".

### **Take care for the use of water based solvents.**

When you are going to use a water based solvent, do not connect to GAS510C and only use the spray dryer.

If need to use the water based solvent when connecting with GAS510C, please refer to P. 36 "About applicable organic solvents".

# 2. Before using this unit

## Precautions when installing the unit



### 1. Always ground this unit

- Be sure to connect the earth wire (the green cable of power cord) to the grounding conductor or ground terminal to prevent electric shock caused by electric leakage.
  - This unit needs a single phase 200-230V power supply. Please entrust the nearest electrical contractor to carry out the works including power connection.
  - The protection impedance of the unit is 0.5Ω or less, ground the unit according to the local technical requirements of the electrical equipment. If the technical requirements are unclear, the grounding works should be constructed and accepted according to the grounding resistance 4Ω or less.
- Do not connect the earth wire to gas or water pipes. If not, fire disaster may be caused.
  - Do not connect the earth wire to the grounding of telephone wire or lightning conductor. If not, fire disaster or electric shock may be caused.

### 2. Please use the special socket. Pay attention to the color of each core wire when connecting the power cord.

<p>Use a socket that match the electric capacity.</p> <p>GAS510C Single phase 200-230V~ 50/60Hz 5.5-12A</p>	Core Wire Color	In-house Wiring
	Black	Voltage Side (L)
	White	Voltage Side (N)
	Green	Ground Side

- Be sure to check that the breaker on the power source equipment side is turned "OFF" when connecting power cord without fail. Note that this unit does not attach the power plug as standard component. Select the appropriate power plug and terminal matching to the power capacity of the power source equipment to be connected, and connect them.
- There could be the case that the unit does not run even after turning ON the power. Inspect whether the voltage of the main power is lowered than the specified value, or whether other device(s) uses the same power line of this unit. If the phenomena might be found, change the power line of this unit to the other power line.
- For connecting of the device to the power source, ask the dealer that you purchased this unit from or an electrical contractor for safe.

### 3. Choose a proper place for installation

Do not install this unit in a place where:

- Flammable gas or corrosive gas is generated.
- Ambient temperature below 5°C or above 35°C.
- Ambient temperature fluctuates violently.
- Place where the water is easy-to-be splashed.

- Rough or dirty surface.
- There is direct sunlight.
- There is excessive humidity and dust.
- There is a constant vibration.



Install this unit on a stable place with the space as shown below.


# Before using this unit

## Precautions when installing the unit


### Warning


#### 4. Do not use this unit in an area where there is flammable or explosive gas

-  Never use this unit in an area where there is flammable or explosive gas. This unit is not explosion-proof. An arc may be generated when the power switch is turned ON or OFF, and fire/explosion may result.
-  Refer to page 59 "15. List of Dangerous Substances".





#### 5. Do not modify

-  Modification of this unit is strictly prohibited. This could cause a failure.




#### 6. Do not topple or tilt this unit

-  Please try to place the unit in a flat place. If the placement is uneven, unexpected accidents may occur.



#### 7. Place the unit

-  Due to sudden earthquake, impact, etc., the product may collapse or move, and then be damaged. It is best to avoid places where there are many people and take safety precautions.

## 2. Before using this unit

### Precautions when installing the unit

#### Warning


#### 8. Pay attention to the use of flammable and explosive solvents.



Take extreme care for use of an explosive or a flammable solvent. Such a solvent may cause an explosion or a fire.

Check whether a solvent can be used or not, please refer to "About the organic solvents" in the section 5. Handling precautions. (P.36)

Always monitor the oxygen concentration in the unit during operation to assure safety.



#### 9. Absolutely prohibit the use of substances that can explode in an oxygen-free environment or substances containing such ingredients



Spray dryer and organic solvent recovery unit are not explosion-proof and can only provide an oxygen-free environment. The use of substances that can explode in an oxygen-free environment or substances containing such ingredients, for example: ammonium nitrate, nitroglycerin, will cause explosion, fire, poisoning or other accidents.

#### 10. The use of substances that release oxygen or contain such ingredients is absolutely prohibited



The use of substances that can release oxygen will destroy the oxygen-free environment and cause accidents such as explosion and fire. For example: peroxides, potassium permanganate.

#### 11. Absolutely prohibit the use of toxic or biohazardous substances




This product and the spray dryer are not developed for biosafety purpose, do not have the ability to treat toxic or biohazardous substances, and the use of toxic or biohazardous substances is absolutely prohibited, for example: polychlorinated biphenyls, cyanide, virus or bacteria.

## 2. Before using this unit



### Precautions when installing the unit

#### Warning

#### 12. Absolutely prohibit the use of substances containing unknown ingredients



-  The thermal decomposition of ingredients with unknown properties may cause explosion, fire, poisoning or other accidents.

#### 13. Handling of power code

-  Do not entangle the power cord. This will cause overheating and possibly a fire.  
Do not bend or twist the power cord, or apply excessive tension to it. This may cause a fire and electrical shock.  
Do not lay the power cord under a desk or chair, and do not allow it to be pinched in order to prevent it from being damaged and to avoid a fire or electrical shock.  
Keep the power cord away from any heating equipment such as a room heater. The cord's insulation may melt and cause a fire or electrical shock.
-  If the power cord becomes damaged (wiring exposed, breakage, etc.), immediately turn off the power at the rear of this unit and shut off the main supply power. Then contact your nearest dealer for replacement of the power cord. Leaving it may cause a fire or electrical shock.  
Connect the power plug to the receptacle which is supplied appropriate power and voltage.



# 3. Names and functions of parts

## Main unit



# 3. Names and functions of parts

## GAS510C Internal mechanism



# 3.Names of parts and their function

## Operation interface overview

### Preheat waiting interface

SINCE 1889  
**Yamato**  
Solvent Recovery Unit  
GAS510C  
Preheating 65 Second Yamato Scientific Co., Ltd. HMI Ver:1.01 PLC Ver:1.01

### Operation curve interface

Graph showing Cooling temp, Oxygen, Blower, and Air volume over time. Values: Cooling temp 0.0°C, Oxygen 0.0%, Blower 0.0%, Air volume 0.00 m³/min.

### Initial interface (install the air flowmeter)

Yamato 2024/04/29 13:13:29  
HMI Ver:1.02 PLC Ver:1.02 DL411C + GAS510C

Blower, Compressor, Small flow N2, Air in, Refrigerator, Defrost, Large flow N2, Exhaust

Blower: PV: 0.0%, SP: 10.0%, Q: 0.09 m³/min  
Cooling temp: PV: 14.6°C, SP: -10.0°C  
Oxygen: PV: 20.5%, SP: 5.0%

Buttons: Run, Alarm, Use nitrogen spray, Manual defrost, Alarm

### Language selection pop-up interface

Yamato 2024/04/28 13:39:43  
HMI Ver:1.02 PLC Ver:1.02 DL411C + GAS510C

Language selection: 中文, 日本語, English

Parameters: PV: 5.0%, SP: 3.0%, Q: 0.09 m³/min

Buttons: Run, Alarm, Use nitrogen spray, Manual defrost, Alarm

### Initial interface (without air flowmeter)

Yamato 2024/04/29 13:11:26  
HMI Ver:1.02 PLC Ver:1.02 DL411C + GAS510C

Blower, Compressor, Small flow N2, Air in, Refrigerator, Defrost, Large flow N2, Exhaust

Blower: PV: 0.0%, SP: 10.0%, No air volume meter  
Cooling temp: PV: 14.3°C, SP: -10.0°C  
Oxygen: PV: 20.6%, SP: 5.0%

Buttons: Run, Alarm, Use nitrogen spray, Manual defrost, Alarm

### System parameters interface

Yamato 2024/04/28 13:39:32  
HMI Ver:1.02 PLC Ver:1.02 DL411C + GAS510C

Oxygen CAL: 0.0%

The system time of GAS is automatically calibrated with the system time of spray dryer.

Heating permissive setpoint: 3.5%, Alarm setpoint: 8.0%

Cooling temp CAL: 0.0°C, SC: 1.000

Refrigerator overload current: In operation: 7.50A, At startup: 10.00A

5 second recording period

Buttons: Save parameters, Alarm

### Alarm interface

Yamato 14:03 Er13. | 2024/04/28 14:05:01  
HMI Ver:1.02 PLC Ver:1.02 DL411C + GAS510C


No.	Alarm date	Time	Alarm information
1	2024/04/28	14:03	Er13. During the use of GAS, the oxygen conce
0	2024/04/28	13:41	Er21. The air volume of GAS is abnormal, and

Buttons: Alarm reset, Home

# 3.Names of parts and their function

## Online operation interface overview (DL411C as an example)


### Preheat waiting interface



Yamato 2024/04/28 14:23:37  
HMI Ver:1.01 PLC Ver:1.01 DL411C + GAS411C GAS411C

Please turn on the power switch of GAS.  
Wait until the GAS warm-up is complete.

### Operation interface




Yamato 2024/04/29 10:26:53  
HMI Ver:1.01 PLC Ver:1.01 DL411C + GAS510C GAS510C

<b>Inlet temp</b> PV: 20.7 °C SP: 180.0 °C Selected Heating Alarm	<b>Outlet temp</b> PV: 24.7 °C SP: 0.0 °C Selected Heating Alarm	<b>Needle</b> PV: 60 sec. SP: 60 sec. Poke twice Manual Auto	<b>Blower</b> MV: 50.0 % Run Alarm
--	---	---	---

Liquid pump: 60.0 rpm

Heater Liquid sending Liquid back Pulse jet Needle Air IN Alarm

### Initial interface




Yamato 2024/04/28 14:17:52  
HMI Ver:1.01 PLC Ver:1.01 DL411C + GAS510C GAS510C

Please assemble the glassware as shown on the right.  
Please insert the outlet temperature sensor into the glassware.  
Please insert the nozzle and connect each pipe.  
After the above is completed, Click the nitrogen in button.  
(If in doubt, please check the instruction manual.)

Nitrogen IN Open the head O2 concentration calibration Alarm

### GAS510C remote control interface




Yamato 2024/04/28 14:20:02  
HMI Ver:1.01 PLC Ver:1.01 DL411C + GAS510C GAS510C

Blower Compressor Small flow N2 Air in  
Refrigerator Defrost Large flow N2 Exhaust

<b>Blower</b> PV: 0.0 % SP: 50.0 % Run Alarm Q: 0.10 m <sup>3</sup> /min	<b>Cooling temp</b> PV: 5.0 °C SP: -10.0 °C Run Alarm Manual defrost	<b>Oxygen</b> PV: 15.9 % Reach Alarm Use nitrogen spray Alarm
--	--	---

### N<sub>2</sub> introducing interface




Yamato 2024/04/28 16:18:22  
HMI Ver:1.01 PLC Ver:1.01 DL411C + GAS510C GAS510C

O<sub>2</sub> : 16.8 % Allowable value: 3.5 %

After the oxygen concentration drops below the allowable value, the Next button will be displayed.  
Click the Next button to enter the running interface.  
Please click the Air IN button to finish.  
After the oxygen concentration returns to more than 18%, the operation can be terminated.

Nitrogen IN Air IN Alarm

### The operation is over and the atmosphere is restored




Yamato 2024/04/29 11:15:41  
HMI Ver:1.01 PLC Ver:1.01 DL411C + GAS510C GAS510C

When the atmosphere is restored, the oxygen concentration must be restored to more than 18% to stop the operation.

End

### Operation jump interface




Yamato 2024/04/29 09:25:36  
HMI Ver:1.01 PLC Ver:1.01 DL411C + GAS510C GAS510C

O<sub>2</sub> : 3.1 % Allowable value: 3.5 %

After the oxygen concentration drops below the allowable value, the Next button will be displayed.  
Click the Next button to enter the running interface.  
Please click the Air IN button to finish.  
After the oxygen concentration returns to more than 18%, the operation can be terminated.

Nitrogen IN Air IN Next Alarm

### Oxygen concentration calibration interface



Yamato 2024/04/28 14:24:50  
HMI Ver:1.01 PLC Ver:1.01 DL411C + GAS510C GAS510C

The oxygen concentration is being corrected for 10 minutes.




9 Minute 28 Second Remaining

# 3.Names of parts and their functions

## Description of switch and indicator lamp in the interface

In the operation interface on the touch screen, the action state of each switch button can be confirmed by checking if the indicator lamp is on. The appearance of the switch button is characterized by a button frame inside which is an effective area for operation.


Type 1: Manual operation button. This kind of button is not only a functional switch but also an indicator lamp for the current state of the function. It has only two states of ON/OFF. The operator will switch the state once every time he operates, and the indicator will switch accordingly. This kind of button has no power-off memory function, and it will be in OFF state after power outage recovery. Its corresponding function can only be activated by the operator.

	Initial state: OFF The temperature controller and heater are not working, and the indicator lamp is in standby state (yellow).
	Click once: the state switches from OFF to ON The temperature controller and heater are working, and the indicator lamp is in operating state (green).
	Click again: the state switches from ON to OFF The temperature controller and heater are not working, and the indicator lamp is in standby state (yellow). The subsequent operations repeat the above actions.

Type 2: Function select button and function enable indicator lamp

This kind of button is only the switch of function, and the usage status of the function is shown by a separate indicator lamp.

This kind of button has power-off memory function, after power outage recovery, will keep the state before power outage. It is characterized by the ability to maintain the last setting of operating parameters.




**Inlet temp**  
PV: 49.8 °C  
SP: 60.0 °C

Selected  
 Heating  
 Alarm

Inlet temp. controller enable button

Status indicator lamp of inlet temp. controller (green means enabled)




**Outlet temp**  
PV: 35.9 °C  
SP: 60.0 °C

Selected  
 Heating  
 Alarm

Outlet temp. controller enable button

Status indicator lamp of inlet temp. controller (green means enabled)



**Needle**  
PV: 33 sec.  
SP: 35 sec.

Poke once  
 Manual  
 Auto

Automatic cleanout needle start/stop button

Status indicator lamp of cleanout needle

Type 3: Status indicator lamp

An indicator lamp indicating the current status of each functional unit. This indicator lamp is internal automatic function and has no corresponding switch.

If the indicator lamp is gray or invisible, it indicates that the function is not running. If the indicator lamp is green or visible, it indicates that the function is running. If the indicator lamp is red, it indicates that the function is in the alarm status.

# 3.Names of parts and their functions

## Value display and input description in the interface

In the operation interface of touch screen, the values can be divided into two types: only display but no input; display and input.

### Type 1: Value display

There is no input box for the value display. Its background and the bottom color are the same, and the display value is white.

The value display will display the data in real time. When an alarm occurs, some specific numerical values will change their colors to prompt the operator.


For example, if the temperature exceeds the upper limit allowed by the equipment, the value will change to red; and the danger level of oxygen concentration is displayed in green, yellow and red colors.

### Type 2: Value input

There is input box for the value input. Its input box is white, and the display value is black.


Click the white input box to pop up the value input popup. Input the required value and press ENTER button to complete the input.

The value input will display the operator's previous input data in real time. Except for the values used in system setting, the inputs of other values are power-off memory type. The operator only needs to operate once in the initial setting, then no operation is required later as long as the value is not changed. (If it is not used for more than 14 days, data may be lost due to PLC internal power supply exhaustion, and may need to be reset.)



- ※ Special note: when the data setting in the system parameter interface is completed, click the Save Parameter button in the lower left, the parameters will be saved only after the button turns green. If the modified data is not saved or the equipment is powered off, the data will be restored to the previous data before modification.

Set the system information correctly at the first startup.



## 4. Operating procedures


### Preparations

(1) Connecting the power cord

First check that the switches of the control assembly and the ELB are OFF and then connect the power cord securely to the power supply meeting the specified voltage and current.

(2) Connecting an earth

The power cord of this unit is earthed 3-core (VCT) including the earth wire, and be sure to ground. The grounding resistance of the standard socket grounding terminal must be less than  $4\Omega$ .




(3) Connecting the exhaust duct

After the operation, the organic solvent, hot air and micro powder in the pipe are discharged through the AIR IN action. Connect the exhaust pipe attached to the exhaust port and exhaust it outside the house using the fume hood. Do not peep into the exhaust port or discharge gas directly into the house, so as to avoid danger. The exhaust port is at the bottom of the left side of the main body.

(4) Connecting the collecting flask

Use a clamp to fix the collecting flask at the connecting port under the ball valve, and then turn on the ball valve after fixing.


If need to remove the collecting flask, turn off the ball valve, then remove the clamp, and finally remove the collecting flask.



## 4. Operating procedures


### Preparations before operation (ADL312SC+GAS510C)

(1) ADL312SC is used on a table. When connecting with GAS, please set the special stand for ADL312SC.



Put the unit body of ADL312SC on the special stand, and use M4×10WS screws to fix the adjusting feet at 4 positions of the fixing plate. Use M4 cross knurled screws to fix the container holder at the front of the special stand.


(2) Remove the suction filter and the blind plate.




# 4. Operating procedures

## Preparations before operation (ADL312SC+GAS510C)

- (3) The attached connecting pipe is installed by M4x10 flat head screws.  
The hose from the cyclone outlet is connected to the installed connecting pipe.



- (4) Connect ADL312SC with GAS510C.




- A: Connect the hose from the heater inlet to GAS510C.  
Put an O-ring and securely fix using a clamp for installation.
- B: Connect the hose from the connecting pipe to GAS510C filter hood inlet.  
Put an O-ring and securely fix using a clamp for installation.
- C: Connect the polyester pipe and fix it with a clamp.
- D: Connect the polyester pipe and fix it with a clamp.
- E: Connect the polyester pipe to the CDA supply source and fix it with a clamp.
- F: Connect the polyester pipe to the N<sub>2</sub> supply source and fix it with a clamp.
- G: Connect an interface cable.
- H: Set GF300 to the main body according to the ADL312SC instruction manual.

## 4. Operating procedures

### Preparation before operation (GB211C+GAS510C)


- (1) Please remove the 4 screws, open the right side door, remove the clamp of the air intake filter, and then remove the air intake filter, as shown in the following figure.




# 4. Operating procedures

## Preparation before operation (GB211C+GAS510C)

- (2) The attached connecting pipe is installed by M4x10 flat head screws.  
The hose from the cyclone outlet is connected to the installed connecting pipe.



- (3) Connect GB211C with GAS510C.




- A: Connect the hose from the connecting pipe to GAS510C filter hood inlet.  
Put an O-ring and securely fix using a clamp for installation.
- B: Connect the hose from the heater inlet to GAS510C.  
Put an O-ring and securely fix using a clamp for installation.
- C: Connect the polyester pipe and fix it with a clamp.
- D: Connect the polyester pipe and fix it with a clamp.
- E: Connect the polyester pipe to the CDA supply source and fix it with a clamp.
- F: Connect the polyester pipe to the N<sub>2</sub> supply source and fix it with a clamp.
- G: Connect the communication port with a network cable.
- H: Set GF301C to the main body according to the GB211C instruction manual.


# 4. Operating procedures

## Preparation before operation (DL411C+GAS510C)


- (1) Remove the fixing screws as shown in the following figure, and remove the online blind plate and the right blind plate of DL411C.



- (2) Remove the clamp and remove the air suction filter




- (3) The attached connecting pipe is inserted in the box from back to front, and then fixed on the panel with the attached double-washer screws M4x10.




## 4. Operating procedures

### Preparation before operation (DL411C+GAS510C)

(4) Fix the GAS510C connecting pipe to DL411C with clamps.



(5) Connect DL411C with GAS510C.



A: Connect the hose from the connecting pipe to GAS510C filter hood inlet.

Put an O-ring and securely fix using a clamp for installation.

B: Connect the hose from the DL411C air inlet to GAS510C.

Put an O-ring and securely fix using a clamp for installation.

C: Connect the polyester pipe and fix it with a clamp.

D: Connect the polyester pipe and fix it with a clamp.

E: Connect the polyester pipe to the CDA supply source and fix it with a clamp.


F: Connect the polyester pipe to the N2 supply source and fix it with a clamp.

G: Connect the communication port with a network cable.

H: Set the drying chamber and glassware to the main body according to the DL411C instruction manual.

# 4. Operating procedures

## Operating method





Turn the ELB on the right side of the main body ON and turn the **Power switch** on the operation panel of the main body ON. Temperature controllers, indicator lamps and the key panel will be displayed.

※The operating method is illustrated by the example of online DL411C and GAS510C. The operation of online ADL312SC and GAS510C, online GB211C and GAS510C are the same.


GAS will take 10 seconds to boot (O<sub>2</sub> concentration meter warm-up time). Please wait patiently to switch to the next interface.

※When GAS does not finish booting, the spray dryer is also in the standby interface.



After GAS is booted, please observe whether the displayed value of O<sub>2</sub> concentration on GAS is above 18% in standby state.

When the O<sub>2</sub> concentration is always lower than 18%, please click the O<sub>2</sub> CAL button on the spray dryer to start the O<sub>2</sub> concentration calibration. After 10-minute calibration, the displayed value of the O<sub>2</sub> concentration should be higher than 18% before proceeding with the following operations.




※When calibrating the O<sub>2</sub> concentration, do not perform other operations and do not power off the equipment. Otherwise, the calibration may fail.

※If the O<sub>2</sub> concentration still cannot exceed 18% after the calibration, please refer to P.49 "The value of oxygen concentration monitor cannot recover beyond 18% after operation".

※ **When the above operations are invalid, considering that the probe of the O<sub>2</sub> concentration meter may have failed, it is forbidden to use this equipment. As for the follow-up treatment, please contact your agent or our customer service.**


# 4. Operating procedures

## Operating method



Please set the required blower output power on the touch screen of GAS. The setting range is 5.0% ~ 100.0%. Please refer to P.30 [ The relation between blower and temperature/dry air amount (reference) ].

※ If need more accurate air volume display, please select an air volume meter to display the real-time air volume changes.



Please set the required condensation temperature on the touch screen of GAS. The setting range is -20.0°C~30.0°C.

The refrigerator contains the delay timing circuit. Because it will take 2mins to start after refrigerator stops, the refrigerator may not activate immediately if press **BLOWER** switch.

**During operation, please according to the characteristics of solvent, set an appropriate control temperature of condenser. Otherwise, when the condenser temperature is too low, the solvent vapor solidifies after liquefaction or directly condenses, resulting in pipeline blockage.**


Please set N<sub>2</sub> supply pressure as 0.1MPa.

After the above preparations are completed, please complete the experimental operation on the spray dryer. The instruction manual takes the operation of DL411C as examples:




(1) Press the **N<sub>2</sub> IN** button to enter the nitrogen lead-in interface.

When the oxygen concentration decreases below the allowable value, the **NEXT** button will appear and click it to enter the running interface.




※ When need the reverse operation to restore the atmosphere, please click the **N<sub>2</sub> IN** button first to turn off the nitrogen lead-in, and then press the **AIR IN** button to jump to the atmosphere recovery interface. Please note that the atmosphere recovery is an irreversible operation. Once entering the atmosphere recovery stage, the operation must be finished before other functions can be used.





# 4. Operating procedures

## Operating method




(2) When the oxygen concentration decreases below the allowable value, the **NEXT** button will appear and click it to enter the running interface.



(3) There are inlet temperature and outlet temperature controllers on the running interface, which can be used for display and temperature setting.

By clicking the icon of inlet temperature or outlet temperature, you can select the inlet temperature or outlet temperature at will. After selecting, the inlet control or outlet control indicator lamp in the inlet temperature or outlet temperature controller will light up.




PV value displays the real-time temp. of the temp. sensor, SP value is black characters on white background, click to set the operating temp.

※ The setting range of each temp. controller is different.

Inlet temp. setting range: 0 ~ 240°C (ADL312SC, GB211C)  
0 ~ 300°C (DL411C)

Outlet temp. setting range: 0 ~ 100°C

e.g.: select the inlet control, set the inlet temp. as 180°C.




(4) Turn ON the Heater switch.


※ After heating, the function switch between inlet temp. control and outlet temp. control will become invalid to prevent the mistake contact in the experiment. If need to switch the control function, firstly turn OFF the Heater switch, switch to inlet temp. control or outlet temp. control, and then turn ON the Heater switch.

# 4. Operating procedures

## Operating method



(5) Set the liquid-sending hose as shown on the left, turn the pull rod of the liquid-sending pump CCW to open the pump head, put the liquid-sending hose in it, and then turn the pull rod CW to make the liquid-sending hose stuck. Insert the other end of the liquid-sending hose into the liquid-sending interface of the spray nozzle. Please use the distilled water as the sample.



ATOMIZING AIR




ON

(6) After the inlet and outlet temp. reach the desired temp., set the spray flow and liquid-sending speed, and turn ON the liquid-sending switch to transport the distilled water.

e.g.: When the outlet temp. reaches about 80°C, the spray flow is set as 10L/min and the liquid-sending speed is set as 60rpm. (**Please refer to the instruction manual of spray dryer.**) Adjust the liquid-sending speed to make the outlet temp. be slightly lower than 75°C.

(7) In order to stabilize the outlet temp. and inlet temp. at the desired temp., please adjust the dry air volume, spray flow and liquid-sending speed again.

e.g.: Adjust the liquid-sending speed to make the outlet temp. be slightly lower than 75°C.



~ Hint ~


- When the inlet temp. is constant, the influences of each setting on the outlet temp. are as follows.  
 Sample liquid-sending volume → small:  
 outlet temp. → high  
 Dry air volume → large: outlet temp. → high  
 Sample concentration (external factor) → high:  
 outlet temp. → high
- If increase the spray flow, the spray droplets will become micronized.
- The volume of spray flow is in direct proportion to the diameter of nozzle orifice.
- When the samples are replaced from the distilled water to the actually used samples, the outlet temp. will become slightly higher due to the non-evaporative part (solid part).

## 4. Operating procedures

### Operating method

(8) When the outlet temp. is stable, replace the samples with the actually used samples. At this point, the outlet temp. will change more or less, if necessary, please adjust the liquid-sending speed again.

e.g.: Replace the samples with 100g sodium chloride 5% solution



(9) During normal spraying, when the sample cannot be sprayed, the orifice of the spray nozzle may be blocked. Click the needle button to squeeze out the blockage.


※ Please refer to the instruction manual of spray dryer.

(10) During normal spraying, if the conical misty samples sprayed from the nozzle becomes irregular, it may be due to the attachment of samples near the orifice of the spray nozzle.

※ Please refer to the instruction manual of spray dryer.

# 4. Operating procedures

## Operating method




~ End process ~



(11) When the sample liquid sending is finished, replace the samples with the distilled water again to clean the nozzle. Clean for about 5mins, turn OFF the liquid-sending switch, and then adjust the spray flow to 0.

e.g.: After about 15mins, when the process of 100g sending liquid is finished, please replace the samples with the distilled water.

(12) Turn the Heater button OFF.



(13) When the inlet temp. is lower than 60°C and the outlet temp. is lower than 50°C, please turn on the **AIR IN** button. Otherwise, you will be prompted:




※ **When the inlet temperature is above 60°C or the outlet temperature is above 50°C, do not stop the operation of the blower by forcibly cutting off the power. This may cause a malfunction, or even lead to combustion or explosion.**

※ At the end of spraying, there is also residual solvent in the pipeline. Please run a little longer to ensure the solvent recovery.

※ When the solvent contains water, there may be frost in the cooling recovery pipe and low solvent recovery. When the solvent recovery in the collecting flask is too low, use the **DEFROST** switch to stop the refrigerator temporarily and defrost the cooling recovery pipe. The refrigerator will automatically be restored after 5 minutes.

(14) Turn off the valve of N<sub>2</sub> supply source when starting to restore the atmosphere.


(15) The END button will appear when the oxygen concentration is restored to more than 18%. Click the END button to restore the standby interface.




## 4. Operating procedures

### Operating method


(16) Please turn OFF the power switch of the spray dryer and GAS.



(17) Remove the container fixing clamp and take out the product collecting container. At this point, please note that the back of the cyclone cover also has powder attached.



(18) Please turn off the hand valve first, then press the collecting flask with your hand, remove the clamp, and finally remove the flask and recover the solvent.




※ The hand valve must be turned off before operation. Because there must be liquid solvent remaining on the inner wall of the condenser and condensing cylinder. If do not turn off the hand valve, the liquid solvent will slowly flow down along the inner wall of the pipe, and finally drip down to the bottom of the box. It may cause flammable gas alarm action, or even will cause a fire.

(19) Wash the containers according to the maintenance method described in the instruction manual of spray dryer.

# 4. Operating procedures

## Gas source selection



The GAS510C internal piston air compressor has the maximum pressure output of 0.3MPa and the maximum flow output of 32L/min. (Note that when there is maximum pressure output, the flow output is 0 L/min, and when there is maximum flow output, the pressure output is 0MPa. The two cannot be reached at the same time.)

If higher pressure or flow are required, use an external gas source to spray.

Operation procedure: on the touch screen of GAS, click the Use nitrogen spray button to adjust the nitrogen gas source to the required pressure.

**※ The maximum pressure of the nitrogen gas source cannot exceed 0.3MPa. Otherwise the pipe may be damaged and leak.**

## 4. Operating procedures

### The relation between blower and temperature/dry air amount (reference)

The following table is the comparison table of each blower scale value and the average dry air amount, please refer to it when operating. If the filter is blocked, the air amount value corresponding to each blower scale value may be reduced. When the value of the filter differential pressure gauge reaches twice the normal value, it needs to be cleaned.

Blower output power (%)	Avg. dry air amount (m <sup>3</sup> /min)
5	0.12
10	0.16
20	0.22
30	0.28
40	0.33
50	0.39
60	0.45
70	0.50
80	0.55
90	0.61
100	0.65

## 4. Operating procedures

### Cooling recovery capacity (reference)

#### 1. Changes in cooling capacity

**!** Since the organic solvent recovery unit uses the air-cooled compressor, it relies on the ambient heat dissipation to provide the cooling capacity, so if the ambient temperature is higher, the cooling capacity will be lower, while if the ambient temperature is lower, the cooling capacity will be higher, the ambient temperature changes, the cooling capacity will also change. **【The higher the ambient temperature, the lower the cooling capacity.】**

The cooling capacity of the organic solvent recovery unit will change with the changes of the ambient temperature and the control temperature. Under a constant condensation temperature condition, when the condensation temperature is lower, the greater the temperature difference with the ambient temperature, the lower the cooling capacity will be.

**【Different condensation temperatures correspond to different cooling capabilities. The lower the condensation temperature, the lower the cooling capacity.】**

#### 2. Formation of cooling recovery capacity

**!** The cooling capacity provided by the cooling unit of solvent recovery unit goes into the closed pipe of cooling recovery, which is divided into three parts:

- ① Reduce the temperature (cooling) of steam and gas in the pipe. The temperature is lowered to the temperature required for steam liquefaction.
- ② To liquefy (recover) the steam. The steam is cooled and liquefied, and the temperature of the gas and liquid is basically unchanged during the liquefaction process.
- ③ Continue to lower the gas and liquid temperature to the set condensation temperature (cooling improves recovery). The cooling capacity provided by the compressor goes into the condensing cylinder through an independent closed pipe. Steam liquefaction occurs on the surface of the cooling pipe. The liquid formed by liquefaction adheres to the surface of the cooling pipe, blocking the heat exchange formed by direct contact between the steam and the pipe, forming a thermal resistance and reducing the efficiency of liquefaction (recovery). At the same time, the gas moves at high speed in the closed pipe, accelerating the gas flow speed on the surface of the liquid substance, and promoting the evaporation of the liquid substance (secondary vaporization). In order to suppress the above adverse effects, it is necessary to set the condensation temperature well below the boiling point temperature of the solvent to achieve the best cooling recovery effect. (The condensation temperature of alcohol solvents is generally lower than 3°C)

#### 3. Effect of inlet temperature of spray dryer on cooling recovery capacity

**!** Assuming that the cooling capacity of the organic solvent recovery unit is certain, other conditions remain unchanged, the higher the inlet temperature, the greater the heat input of the spray dryer, the greater the temperature difference formed with the set condensation temperature, the more the cooling capacity required for the cooling process, the refrigeration capacity left for cooling to improve the recovery rate will be reduced, and the actual condensation temperature will increase.

The organic solvent recovery unit has the automatic adjustment function of cooling capacity. After the inlet temperature is adjusted to be higher, the organic solvent recovery unit will automatically adjust the cooling capacity to keep the condensation temperature constant. When the input heat exceeds the limit cooling capacity of the organic solvent recovery unit, **【the actual condensation temperature begins to rise slowly from the set condensation temperature, at this time the recovery rate will decrease】**.

## 4. Operating procedures

### Cooling recovery capacity (reference)

#### 4. Effect of liquid sending amount on cooling recovery capacity

**!** Under the same conditions, the greater the liquid sending amount, the more the solvents need to be condensed in the same time, and the more the needed cooling capacity. Once exceeding the limit cooling capacity of the organic solvent recovery unit, [the actual condensation temperature begins to rise slowly from the set condensation temperature, at this time the recovery rate will decrease].

#### 5. Effect of circulating gas volume on cooling recovery capacity

**!** The effect of circulating gas volume on cooling recovery capacity is mainly reflected in the following three aspects:

- ① Under the same conditions, the greater the gas volume, the more the hot gas flowing per unit time, to make more hot gas cooling requires more cooling capacity, the refrigeration capacity left for cooling to improve the recovery rate will be reduced.
- ② The greater the gas volume, the faster the steam flows through the condensing cylinder, the shorter the time, greatly reducing the opportunity and time of contact between the steam and the cooling pipe, the liquefaction is not sufficient, and the recovery rate will be reduced.
- ③ The greater the gas volume, the faster the gas flow on the surface of the liquid substance, which promotes the evaporation of the liquid substance (secondary vaporization).

In summary, [the greater the gas volume, the weaker the cooling recovery capacity]. This is also the reason why the maximum gas volume of the organic solvent recovery unit is lower than the maximum gas volume of the standalone spray dryer.

#### 6. Effect of solvent type on cooling recovery capacity


**!** The solvent types are different (including different ratios), then the solvent density, viscosity, specific heat, boiling point, freezing point, latent heat of vaporization, vaporization volume, thermal conductivity are different, need to be treated differently.

Take anhydrous ethanol and water as an example: when the liquid sending speed is the same, ideally, 1L anhydrous ethanol and 1L water are input to the spray drying and cooling recovery system under the same conditions. When the condensation temperature is 0°C, the anhydrous ethanol can be recycled well. The condensation temperature of water cannot be lower than 0° C, otherwise the water will freeze on the cooling pipe, the recycle cannot be performed properly.

# 4. Operating procedures


## Cooling recovery capacity (reference)

The following chart shows the corresponding curve between the output power of Roots blower and the maximum operating temperature of the inlet, under the condition of a certain liquid sending amount and based on ethanol as the standard liquid. (When the ambient temperature is 20°C, the power supply voltage is 220V, the spray is continued for 1 hour, and the actual condensation temperature of the organic solvent recovery unit is maintained below 3°C) Please use it as reference during operation. [The area surrounded by curves and coordinate axes is the usable area]



# 4. Operating procedures


## Cooling recovery capacity (reference)




# 5. Handling Precautions

## Warning

### 1. Substances that cannot be used


-  Such substances may cause an explosion or a fire. Whether a solvent may be used or not shall be judged according to "About specified organic solvents" in section 5. Handling precautions.  
Always monitor the oxygen concentration during operation to assure safety.  
See "15. List of Dangerous Substances" on P.59.

### 2. If a problem occurs


-  If smoke or strange odor should come out of this unit for some reason, turn off the power key right away, and then turn off the circuit breaker and the main power. Immediately contact a service technician for inspection. If this procedure is not followed, fire or electrical shock may result. Never perform repair work yourself, since it is dangerous and not recommended.

## Caution


### 1 Do not put anything on this unit.

-  Do not put anything on this unit. It will cause injury if fall.


### 2. During a thunder storm

-  During a thunderstorm, turn off the power key immediately, then turn off the circuit breaker and the main power. If this procedure is not followed, fire or electrical shock may be caused.



### 3. After installing

-  It may cause injury to a person if this unit falls down or moves by the earthquake and the impact, etc.. To prevent, take measures that the unit cannot fall down.

### 4. Clean the collecting flask during operation.

-  When the solvent collected in the collecting flask is about to be filled, because the condensing cylinder has a certain capacity, the collecting flask can be cleaned without stopping operation.  
Please turn off the hand valve at the connecting pipe of the collecting flask before removing the collecting flask. Safely dispose the solvent in the collecting flask as soon as possible, then reinstall the collecting flask and turn on the hand valve. The whole operation should be controlled within 5 minutes.  
If the solvent is a special reagent, the transfer is hard and the disposal time is long, please order another collecting flask from our company or agent to replace it for safety.

### 5. It is strictly prohibited to remove the collecting container to collect samples during operation.

-  It is strictly prohibited to remove the collecting container during operation which will lead to the rapid inhalation of a large amount of air in the sealed pipeline, and there is high temperature in the pipeline, the residual solvent vapor is easy to cause explosion when encountering the oxygen.
-  Samples should be collected after the end of the operation.

# 5. Handling precautions

## About applicable organic solvents

This unit has been designed to use isopropyl alcohol and ethanol. Note that the following restrictions shall apply when other organic solvents are used.

(1) Restrictions because of explosion limit oxygen concentration

When N<sub>2</sub> gas, organic solvents, and air (oxygen) are mixed, explosion will occur when the oxygen concentration is over the oxygen concentration limit and if there is an ignition source. This means that it is desirable the oxygen concentration limit of a solvent is as high as possible.

The oxygen concentration limit for this unit is 9%. Do not use any organic solvent whose oxygen concentration limit is below 9%. (See "Oxygen concentration limit table" on P.25. Also, calculate limits for organic solvents not shown in the table using the "calculation method of concentration limits".

(2) Restrictions because of the boiling point

Although the unit recovers solvent by cooling it in the condenser (cooling trap), the outlet temperature of the condenser (TRAP TEMP MONITOR) may become considerably high depending on the outlet temperature, amount of liquid that is sent, environmental temperature, or duration of sending liquid.

[Example]

Environmental temperature: 35°C	Amount of liquid to send: 2000mL/H (ethanol)
Air amount: 0.5m <sup>3</sup> /min	Liquid sending time: 30 min
Outlet temperature: 100°C	Condenser outlet temperature: 37°C

Thus a solvent with a lower boiling point may not be condensed efficiently. (See "Table of oxygen concentration limit" on P.25.) In fact, you do not need to set inlet and outlet temperature higher when you use a solvent with a low boiling point. For example, you can operate the unit with lower condenser outlet temperature by reducing the air amount and amount of liquid to send for solvents with a low boiling point such as methylene chloride.

[Example]

Environmental temperature: 25°C	Amount of liquid to send: 1170mL/H (methylene chloride)
Air amount: 0.45m <sup>3</sup> /min	Liquid sending time: 20 min
Outlet temperature: 38°C	Condenser outlet temperature: 14°C

(3) Restrictions because of boiling points

When the solvent contains an organic solvent with a high melting point or water, too low trap temperature might cause it condensate in the condenser. In such a case, raise the trap temperature by increasing air amount, increasing the inlet temperature, or increasing the amount of liquid to send, or stop the freezer once with DEFROST switch and allow condensed solvent or water to solve. (When the trap temperature is lower than the melting point before sending solvent, watch the trap temperature a while because that temperature will rise once sending of solvent is started.)

When you use a water soluble solvent, try to disconnect this unit and operate with ADL312SC or GB211A only.

(4) Restrictions because of corrosion resistance

This unit has been designed to use isopropyl alcohol and ethanol. When other solvents are used, care must be taken because service lives will differ from part to part. See the table of corrosion resistance on P.26 and the piping system drawing that shows parts made of materials other than stainless steel or glass on P.27. In terms of this unit, components other than the liquid sending tube are exposed to thin solvent steam only and will not immediately corrode even if they are exposed to chemicals marked with  $\Delta$  or x in the corrosion resistance table, if any abnormalities such as abnormal increase speed of oxygen concentration or a gas leakage in the pipe path, replace defective parts immediately because service life of some parts may be shortened due to solvents other than ethanol and isopropyl alcohol.

## 5. Handling precautions

### About applicable organic solvents

Chemicals that each type of liquid sending tubes are as follows.

Silicone tube: ethanol, isopropyl alcohol, methanol, acetone, acetic ether

Viton tube: xylene, toluene, benzene, hexane, chloroform, methylene chloride

[Oxygen concentration limit table]

Organic solvent	Boiling point [°C]	Melting point [°C]	Oxygen concentration limit [%]
Xylene	(o) 144 (m) 139.3 (p) 138.5	(o) -25 (m) -47.4 (p) 13.2	(o) 10.5 (m) 11.5 (p) 11.5
Isopropyl alcohol	82.3	-88	9.0
Benzene	80.1	5.5	10.5
Ethanol	78.4	-114.3	9.9
Acetic ether	77.1	-83.6	10.0
Hexane	67.7	-95.3	11.4
Methanol	64.6	-97.4	9.7
Chloroform	61.2	-63.5	Non-combustible
Acetone	56.2	-94.6	10.4
Methylene chloride	40	-97.7	23.9

[How to calculate an oxygen concentration limit]

Molecular formula of flammable gases:  $C_a H_b O_c N_d S_e F_f$  (F means halogen) d is not used.

$$v = \frac{100}{1 + 4.773 \left[ a + e + \frac{b-2c-f}{4} \right]}$$

Oxygen concentration limit

$$= [100 - \{L + (1 - L/v) \times 100\}] \times 0.209 \quad [\%]$$

L: Lower explosion limit of a flammable gas [%]

Example: Toluene

According to the molecular formula  $C_6H_5CH_3$ ,  $a=7$ ,  $b=8$ ,  $c=e=f=0$

Lower explosion limit  $L=1.2$

$$v = \frac{100}{1 + 4.773 \left[ 7 + 0 + \frac{8-0-0}{4} \right]} = 2.27$$

Oxygen concentration limit

$$= [100 - \{1.2 + (1 - 1.2/2.27) \times 100\}] \times 0.209$$

$$= 10.7\%$$

# 5. Handling precautions

## Corrosion resistance table

○: Usable    △: Avoid using preferably    x: Unusable for use

Material	Silicone rubber	Viton (FPM)	Chloroprene rubber (CR)	Nitrile rubber (NBR)	Steel acrylic phthalic acid resin paint	POM	Phenol	Polypropylene (PP)	Vinyl chloride (hard)	Polyacetal (PA)
Applicable parts	Glass connecting packing Liquid sending tube Diaphragm cap for the differential pressure meter Bond to glass	O-ring Solenoid valve seal Liquid sending tube Hoses Packing for installing the oxygen concentration measurement sensor	Filter bottom packing Oxygen concentration meter (pump valve, diaphragm)	Nozzle packing Blower oil seal Pressure meter packing Needle valve packing	Blower air contact part	Needle valve BOX for installing the oxygen concentration measuring sensor	Bonding of aluminum honeycomb Hoses	Compressor cover		Tube coupler
Xylene	x	○	x	△~x	△	○	○	△	x	○
Toluene	x	○	x	△~x	△	○	○	△	x	○
Isopropyl alcohol	○	○	○	○	○	△	○	○	○	○
Benzene	△	○	x	x	△	△	○	△	x	○
Ethanol	○	○	○	○	○	○	○	○	○	○
Acetic ether	△	○	△	△~x	△	○	○	△	x	⊗
Hexane	x	○	○	○	△	○	○	△	○	⊗
Methanol	○	x	○	○	○	△	△	○	○	△
Chloroform	x	○	x	x	△	x	○	x	x	x
Acetone	○~△	x	○~△	△~x	△	△	○	△	x	○~△
Methylene chloride	x	○	x	x	△	x	○	△	x	x
Ethylene chloride	○~△	○	x	x	△	x	⊗	○	x	x

# 5. Handling Precautions

## Precautions during operation

- (1) Be sure to connect the earth wire when connecting power supply.
- (2) Securely connect with the spray dryer.
- (3) Using the blower at a low air amount may cause a malfunction. Set it at least 0.2m<sup>3</sup>/min or more.
- (4) Flow meter for the oxygen sensor (small) and a flow meter for introducing N<sub>2</sub> (large) are located inside the left side door. Adjust the flow meter for the oxygen sensor to 0.2L/min at AIR IN during preparations and adjust the flow meter for introducing N<sub>2</sub> to 30L/min when performing N<sub>2</sub> IN.
- (5) According to the connection with the spray dryer, the blower of spray dryer is not used. Please set the air amount on GAS.
- (6) Make sure that the glass chamber is fixed at the specified position without any play. Check the connecting assembly when oxygen concentration will not decrease properly or its increase speed has risen. When a defect at a connection of pipes or at a packing is suspected, repair or replace the relevant part.
- (7) The outlet temperature should not exceed 140°C, because the material of the suction/exhaust hose, material of the filter and the performance of the blower may be deteriorated..
- (8) Check the glass chambers are fixed to the specified position with no gap, and then turn on the switches of blower and heater.
- (9) The unit is not explosion proof. You cannot use this unit in an atmosphere where a flammable gas exists or at a place where a flammable gas may be produced.
- (10) If there is a gap between the product collecting container and the bracket at lower of the cyclone, the dried powder may be accumulated at the lower of the cyclone without falling into the product container. Therefore, pay special attention when installing product container.
- (11) When the heater is ON, do supply the air to the heater part for at least 0.2m<sup>3</sup>/min.
- (12) If the sample is not fed from the feeding pump, the following causes may be considered; the sample tube is crushed at the roller of the pump, the inner wall of the tube is adhered tightly without restoration, or the inner of the nozzle is blocked. Remove the cause, and reset to the normal status.
- (13) Do not perform unattended operation during activating the unit. Since the unit is in idling status and the nozzle is blocked of after the operation using sample, the temperature around outlet is increased and the remaining sample is flown from the sample tube disconnected from the unit, and these failures may cause the indeterminism accident.
- (14) There are two types of specimen tubes, those made of silicone and those made of Viton. Take care they might be corroded, swell and break with some solvents during operation.

## 5. Handling Precaution

### Caution during operation

---

- (15) When the high temperature is set to the temperature around inlet for the operation, supplying too excessive airflow of the blower to the unit may not reach the temperature to the setting one caused by not keeping balance with the heater capacity. To resolve this error, turn down the airflow of the blower, increase the setting temperature, and operate this unit. To avoid a malfunction of the blower, set the air amount below the red zone on the dial scale. When you operate the unit at a higher set temperature, the setting and the actual inlet temperature will not match.  
Moreover, the heater will automatically stop when the inlet temperature exceeds 260°C (320°C for DL411C) and the outlet temperature exceeds 140°C.
- (16) If this unit is not operated, turn "OFF" the earth leakage breaker on the back of the unit.
- (17) The cyclone may charge easily with static electricity depending on the specific specimen used, or operating environment or conditions. Implement countermeasures against static electricity such as attaching included earth clips at three positions on the clamp at the connection of the cyclone or attaching an antistatic brush to the body of the cyclone.
- (18) Depending on the sample to be processed, the static electricity may be occurred at cyclone. Therefore, remove the static electricity with an appropriate method. It is efficient that the wire is wound to the glass portion for grounding, but it is more convenient to use the static electricity remover by setting against the cyclone vertically.
- (19) The filter for the oxygen sensor is consumable. Please replace it at least once every six months. The filter is inside the left door of the equipment.

# 6. Maintenance Method

## Daily Inspection and Maintenance

### Warning

- Disconnect the power cable from the power source when doing an inspection or maintenance unless needed.
- Perform the daily inspection and maintenance after returning the temperature of this unit to the normal one.
- Do not disassemble this unit.


### Caution

- Use a well-drained soft cloth to wipe dirt on this unit. Do not use benzene, thinner or cleanser for wiping. Do not scrub this unit. Deformation, deterioration or color change may result in.



### Monthly maintenance

- Check the earth leakage breaker function.
  - Connect the power cord.
  - Turn the breaker on.
  - Push the red test switch by a ballpoint pen etc. If there is no problem, the earth leakage breaker will be turned off.



# 6. Maintenance Method

## Daily Inspection and Maintenance

 **Caution**


### Filter cleaning

- Filter element

The clogging standard of filter mesh is that the indication on the filter differential pressure gauge is about twice of ordinary indication.

When the filter is clogged, open the left side door of the unit, open the filter case cover, tear off the filter seal sticker, and then loosen and remove the butterfly nut that fix the filter. When assembling, please paste the accessorial filter seal sticker in the last step, after pasting, the concave of filter will be completely hidden.

To remove clogging of filter mesh, beat lightly to remove dusts off. (The filter cannot be washed by water.) Replace with a new one if clogging occurs too often after cleaning.



- Filter case

Remove the filter element and remove dusts attached inside the filter case or on the bottom by suctioning with an electric cleaner or by wiping out with a cloth. During this work, take care so that dusts will not fall on the pipe port on the bottom.

- Refrigerator filter

There is a filter for the freezer capacitor at the lower right position of the front of the unit. The filter cover is held with magnets. Lightly pull it toward you by hand to remove it and then clean the filter. To clean, gently wash it with water or remove dusts using a cleaner. When capacitor fins at the deep back of the filter are clogged with dusts, remove them with a cleaner.

- Filter (for oxygen)

Open the door at the left side of the unit and remove the filter fixed to the front pole.

You can remove the filter holding screws by hand. (See P.9 of the operation manual.)

Rough replacement interval of the filter shall be about six months.

# 6. Maintenance procedures

## Daily inspection / care

### Supply oil to the blower


Exchange lubrication oil after 500 hours of operation (about 20 days) of initial operation after installation and each of about 2000 hours of operation (about three months) thereafter.

- ① Remove the blind plates on the right side plate and at the lower right of the rear side.
- ② Stop the blower, remove the oil drain plug, and allow all old oil to discharge.
- ③ Replace the plug, and fill new lubrication oil to the line of the oil gauge taking care dusts will not enter at the upper oil filler port. (Approx.0.3L)

#### Lubrication oil brands

Appropriate oil is additional turbine oil No.3 (#180) of JIS K2213 or equivalent. Typical brands are as follows.

	Mobil	Esso	Idemitsu	Mitsubishi	Nisseki	Showa Shell	Kyoseki	Cosmo
No. 3	DTE oil heavy medium	Teresso 68	Daphne turbine oil 68	Diamond turbine oil 68	FBK turbine 68	Turbo oil T68 J-H oil68	Kyoseki RIX turbine 68	Cosmo turbine super 68



## 7. Long storage and disposal

### When not using this unit for long term / When disposing

#### Caution

##### When not using this unit for long term...

- Turn off the earth leakage breaker and original power source for safe without fail. Also, store the glass unit after removing it from the main unit. When the glass unit is contacted to the external, it may cause the breakage.

#### Warning

##### When disposing...

- Keep out of reach of children.
- Remove the power cord.

### Matters to consider when disposing of the unit

Environmental protection should be considered

- We request you to disassemble this unit as possible and recycle the reusable parts considering to the environmental protection. The feature components of this unit and materials used are listed below.

Component Name	Material
<b>Parts of Main Unit</b>	
Exterior	Bonderizing steel plate baked with melamine resin coating, stainless steel
Insulating material	PE foam board, ceramic fiber cotton
Condenser (evaporator) Filter case Pipes and joints	Stainless steel
Label	Polyethylene (PET) resin film
Hose	Silicon rubber, Teflon, Viton
<b>Electrical Parts</b>	
Refrigerator	Stainless steel, iron, copper, aluminum, etc.
Compressor	Iron, PP, etc.
Circuit boards	Composites with board, condenser, resister and transformer
Power cord & wiring materials and others	Synthetic rubber, resins
Sensor	Stainless steel and others

# 8. When a trouble occurs

## Safety unit and error indications

The table shows possible causes of activation of the safety unit and solutions.

[Error indication]

When an abnormality occurs to the inlet temperature controller or the outlet temperature controller, the touch screen at the operation panel displays the error screen. When an abnormality occurs, confirm the error content and implement appropriate solutions.

Display	Reasons	Solutions
Er01. The PLC analog module of GAS is faulty.	<ul style="list-style-type: none"> <li>① The wire connection of the PLC analog module is loose</li> <li>② The PLC analog module is damaged</li> </ul>	<ul style="list-style-type: none"> <li>① Power off and restart</li> <li>② If it cannot reset after power off and restart, please contact our service department or agent.</li> </ul>
Er02. GAS temperature transmitter disconnection alarm.	<ul style="list-style-type: none"> <li>① The wire connection of cooling temperature transmitter is loose</li> <li>② The cooling temperature transmitter is damaged</li> </ul>	<ul style="list-style-type: none"> <li>① Power off and restart</li> <li>② If it cannot reset after power off and restart, please contact our service department or agent.</li> </ul>
Er03. The cooling temperature sensor of the GAS is malfunctioning.	<ul style="list-style-type: none"> <li>① The wire connection of cooling temperature sensor is loose</li> <li>② The cooling temperature sensor is damaged</li> </ul>	<ul style="list-style-type: none"> <li>① Power off and restart</li> <li>② If it cannot reset after power off and restart, please contact our service department or agent.</li> </ul>
Er04. The cooling temperature controller of GAS alarms.	<ul style="list-style-type: none"> <li>① The cooling temperature controller data is abnormal</li> <li>② The PLC internal data is confused or lost</li> </ul>	<ul style="list-style-type: none"> <li>① Power off and restart</li> <li>② If it cannot reset after power off and restart, please contact our service department or agent.</li> </ul>
Er05. The Roots blower of GAS alarms.	<ul style="list-style-type: none"> <li>① The current of Roots blower is too large</li> <li>② The Roots blower is damaged</li> </ul>	<ul style="list-style-type: none"> <li>① Please set the blower power to 80% and restart after power off (the current is large when the blower starts at low speed).</li> <li>② If it cannot reset after power off and restart, please contact our service department or agent.</li> </ul>
Er06. The current transmitter of the refrigeration unit of GAS is disconnected.	<ul style="list-style-type: none"> <li>① The wire connection of refrigerating unit current transmitter is loose</li> <li>② The refrigerating unit current transmitter is damaged</li> </ul>	<ul style="list-style-type: none"> <li>① Power off and restart</li> <li>② If it cannot reset after power off and restart, please contact our service department or agent.</li> </ul>

## 8. When a trouble occurs

### Safety unit and error indications

Display	Reasons	Solutions
Er07. Refrigerator overload of GAS.	<ul style="list-style-type: none"> <li>① The air inlet and outlet of compressor are blocked</li> <li>② The condensation temperature is maintained at a relatively high temperature for a long time.</li> <li>③ The compressor is damaged</li> </ul>	<ul style="list-style-type: none"> <li>① Clean the air inlet and outlet</li> <li>② Reduce the liquid sending amount or reduce the air speed</li> <li>③ If it cannot reset, please contact our service department or agent.</li> </ul>
Er08. The operation current of the refrigeration unit of GAS is too small to operate normally.	<ul style="list-style-type: none"> <li>① The compressor starting current is large, the thermal overload protection cuts off the starting circuit.</li> <li>② The compressor is damaged</li> </ul>	<ul style="list-style-type: none"> <li>① Power off and restart</li> <li>② If it cannot reset after power off and restart, please contact our service department or agent.</li> </ul>
Er09. A flammable gas leak alarm occurred in GAS.	<ul style="list-style-type: none"> <li>① A flammable gas leakage is detected</li> <li>② The flammable gas detector is damaged</li> </ul>	<ul style="list-style-type: none"> <li>① Immediately stop operation, power off, and ventilate the experimental area.</li> <li>② If it cannot reset, please contact our service department or agent.</li> </ul>
Er10. The evaporator temperature of GAS is abnormal. (Over 40°C)	<ul style="list-style-type: none"> <li>① The condensation temperature exceeds 40°C</li> <li>② The temperature sensor is damaged</li> </ul>	<ul style="list-style-type: none"> <li>① Reduce temperature, air speed and liquid sending amount</li> <li>② If the evaporator cannot be reset after the temperature is reduced, please contact our service department or agent.</li> </ul>
Er11. The oxygen concentration meter of GAS is disconnected.	<ul style="list-style-type: none"> <li>① The wire connection of oxygen concentration meter is loose</li> <li>② The oxygen concentration meter is damaged</li> </ul>	<ul style="list-style-type: none"> <li>① Power off and restart</li> <li>② If it cannot reset after power off and restart, please contact our service department or agent.</li> </ul>
Er12. The oxygen concentration meter of GAS is abnormal. Oxygen concentration is lower than 0% or higher than 25%.	<ul style="list-style-type: none"> <li>① The oxygen concentration is below 0% or above 25%</li> <li>② The service life of oxygen concentration probe is exhausted</li> </ul>	<ul style="list-style-type: none"> <li>① Restart after power off</li> <li>② If it cannot reset, please contact our service department or agent.</li> </ul>

## 8. When a trouble occurs

### Safety unit and error indications

Display	Reasons	Solutions
Er13. During the use of GAS, the oxygen concentration rose abnormally.	<ul style="list-style-type: none"> <li>① N<sub>2</sub> is exhausted</li> <li>② There are gaps between glassware and pipes</li> <li>③ Remove the collecting flask without turning off the hand valve of collecting flask</li> <li>④ The product collecting container is removed</li> </ul>	<ul style="list-style-type: none"> <li>① Replace with a new gas source</li> <li>② Power off and reconnect the glassware and pipes</li> <li>③ Please turn off the valve before removing the collecting flask</li> <li>④ It is strictly prohibited to collect samples during operation</li> <li>※ If the spray dryer is used separately, the continuous sample recovery operation can be performed, but once GAS510C is used, such operation must be prohibited, because once the air is mixed into the pipeline, it may cause explosion.</li> </ul>

When select the optional air flowmeter, the following table shows the reasons and solutions when the safety device activates:

Display	Reasons	Solutions
Er20. The air volume meter of GAS is disconnected.	<ul style="list-style-type: none"> <li>① The wire connection of air flowmeter is loose</li> <li>② The air flowmeter is damaged</li> </ul>	<ul style="list-style-type: none"> <li>① Restart after power off</li> <li>② If it cannot reset, please contact our service department or agent.</li> </ul>
Er21. The air volume of GAS is abnormal, and the air volume is too small during operation.	<ul style="list-style-type: none"> <li>① The inlet and outlet of air flowmeter are inverted</li> <li>② The air flowmeter is installed at the GAS inlet pipe</li> <li>③ The pipe or filter is clogged</li> <li>④ The air flowmeter is damaged</li> </ul>	<ul style="list-style-type: none"> <li>① Please install in the direction of air flow</li> <li>② Please install at the GAS outlet pipe</li> <li>③ Clean the pipe or filter</li> <li>④ After power off and restart, if it cannot reset, please contact our service department or agent.</li> </ul>

# 8. In the Event of Failure...

## Trouble Shooting

Symptoms	Possible causes	Countermeasures
The POWER does not turn ON	<ul style="list-style-type: none"> <li>● ELB is turned OFF</li> <li>● Malfunction of the power supply</li> <li>● The power cord is disconnected</li> <li>● Malfunction of power switch</li> </ul>	<ul style="list-style-type: none"> <li>● Turn the ELB ON</li> <li>● Check the power supply circuit</li> <li>● Replace the power cord</li> <li>● Replace the power switch</li> </ul>
Although GAS is connected, the operation interface of touch screen does not switch	<ul style="list-style-type: none"> <li>● Defective wires connection</li> <li>● Defective wires</li> </ul>	<ul style="list-style-type: none"> <li>● Connect correctly as per the instruction manual</li> <li>● Replace the wires</li> </ul>
Blower does not activate.	<ul style="list-style-type: none"> <li>● Incorrect connecting of the connector of blower</li> <li>● Lubricating oil level is low</li> <li>● Breaking of blower input cord</li> <li>● Blower switch failure</li> <li>● Blower or inverter failure</li> <li>● Blower circuit failure and wiring failure</li> </ul>	<ul style="list-style-type: none"> <li>● Connect correctly.</li> <li>● Exchange oil periodically</li> <li>● Replace the cord.</li> <li>● Replace the touch screen, PLC or temperature controller.</li> <li>● Replace the blower or inverter</li> <li>● Maintain or replace the part</li> </ul>
The refrigerator does not operate.	<ul style="list-style-type: none"> <li>● The delay timer is in operation</li> <li>● Disconnection of the heater</li> <li>● Touch screen or PLC is defective</li> <li>● Refrigerator circuit and wiring are defective</li> </ul>	<ul style="list-style-type: none"> <li>● Wait for 5 minutes and check whether the refrigerator operates or not</li> <li>● Replace the heater</li> <li>● Replace the touch screen or PLC</li> <li>● Repair the defective parts or replace the temperature controller</li> </ul>
Spray pressure will not increase.	<ul style="list-style-type: none"> <li>● Dial is set at 0</li> <li>● Touch screen or PLC is defective</li> <li>● Imperfect connection of the hose</li> <li>● Compressor is defective</li> <li>● Defective solenoid valve</li> <li>● Compressor circuit or wiring is defective</li> </ul>	<ul style="list-style-type: none"> <li>● Adjust the dial</li> <li>● Replace the touch screen or PLC</li> <li>● Repair or replace the defective points</li> <li>● Replace</li> <li>● Replace</li> <li>● Repair the defective point</li> </ul>
Oxygen concentration will not decrease.	<ul style="list-style-type: none"> <li>● N2 supply source is defective</li> <li>● Flow meter for introducing N2 is clogged</li> <li>● Imperfect connection of the hose</li> <li>● Defective solenoid valve</li> <li>● Touch screen or PLC is defective</li> </ul>	<ul style="list-style-type: none"> <li>● Adjust appropriately</li> <li>● Replace the flow meter</li> <li>● Repair or replace the defective points</li> <li>● Replace the solenoid valve</li> <li>● Replace the touch screen or PLC</li> </ul>
The value of oxygen concentration monitor cannot go beyond 18%	<ul style="list-style-type: none"> <li>● High altitude, low oxygen content</li> <li>● Flowmeter of oxygen sensor is blocked</li> <li>● Pump failure</li> <li>● Filter degraded or blocking</li> <li>● Oxygen sensor failure</li> <li>● Oxygen sensor board failure</li> </ul>	<ul style="list-style-type: none"> <li>● Use it at proper altitude or use default oxygen concentration</li> <li>● Replace the flowmeter</li> <li>● Replace</li> <li>● Replace the filter</li> <li>● Replace</li> <li>● Replace</li> </ul>
The value of oxygen concentration monitor does not change	<ul style="list-style-type: none"> <li>● Flowmeter of oxygen sensor is blocked</li> <li>● Pump failure</li> <li>● Filter degraded or blocking</li> <li>● Oxygen sensor failure</li> <li>● Oxygen sensor board failure</li> </ul>	<ul style="list-style-type: none"> <li>● Replace the flowmeter</li> <li>● Replace</li> <li>● Replace the filter</li> <li>● Replace</li> <li>● Replace</li> </ul>

# 8. In the Event of Failure...

## Trouble Shooting



Symptoms	Possible causes	Countermeasures
The value of oxygen concentration monitor cannot recover beyond 18% after operation	<ul style="list-style-type: none"> <li>● Influenced by active carbon filter, the internal air oxygen concentration cannot recover to atmosphere oxygen concentration</li> <li>● Flowmeter of oxygen sensor is blocked</li> <li>● Pump failure</li> <li>● Filter degraded or blocking</li> <li>● Oxygen sensor failure</li> <li>● Oxygen sensor board failure</li> </ul>	<ul style="list-style-type: none"> <li>● Use gas pipe switching valve, connect the test port of oxygen concentration meter to atmosphere directly. See pictures on this page.</li> <li>● Replace the flowmeter</li> <li>● Replace</li> <li>● Replace the filter</li> <li>● Replace</li> <li>● Replace</li> </ul>
Temperature controller failure	<ul style="list-style-type: none"> <li>● Defective display function</li> <li>● Sensor failure</li> </ul>	<ul style="list-style-type: none"> <li>● Maintain the defective points or replace the PLC</li> <li>● Replace the sensor</li> </ul>
Failure of blower power regulation	<ul style="list-style-type: none"> <li>● Adjusting circuit failure and wiring failure</li> <li>● Insufficient heater capacity due to excessive dry air flow</li> <li>● Invalid blower power setting for spray dryer or GAS510C</li> </ul>	<ul style="list-style-type: none"> <li>● Maintain the defective points or replace the PLC</li> <li>● No error. For operating this unit with high temperature, decrease the flow rate of the dry air or increase the setting value</li> <li>● When the blower power setting of the spray dryer is invalid, try to use the blower power setting of GAS510C. The converse is equally valid.</li> </ul>


◆ In case if the error other than listed above occur, please immediately cut off the power switch of the unit body, unplug the power cord, and contact the sales agent, our company's business office or customer service center.

※ When the value of oxygen concentration monitor cannot recover beyond 18% after operation, it's suspected that it may be influenced by active carbon filter adsorption, operate as per the following pictures. Adjust the gas pipe switching valve to connect to atmosphere (External gas), wait for several minutes, and check if the oxygen concentration display recovers beyond 18%. After confirmation, make reverse operation to recover the working state (Internal gas).

Open the access door at left side



Turn the gas pipe switching valve from  to 



## 9. After Service and Warranty

### When requesting a repair

#### When requesting a repair

If any trouble occurs, immediately stop operation, turn the power switch off, pull out the power plug and contact your dealer, our sales office or our customer service center.

Information necessary for requesting a repair

- Model name of the product
- Serial number
- Date (y/m/d) of purchase
- Description of trouble  
(as in detail as possible)

} See the warranty card or the nameplate on the unit.  
See the section "3. Names and functions of parts" on page 9.

Be sure to indicate the warranty card to our service representative.

#### Warranty card (attached separately)

- Warranty card is given by your dealer or one of our sales offices and please fill in your dealer, date of purchase and other information and send it to our customer service center by Facsimile (03-3231-6523). Then, store it securely.
- Warranty period is one full year from the date of purchase. Repair service for free is available according to the conditions written on the warranty card.
- For repairs after the warranty period consult your dealer, one of our sales offices or our customer service center.  
Paid repair service is available on your request when the product's functionality can be maintained by repair.

#### Minimum holding period of repair parts

The minimum holding period of repair parts for this product is seven years after end of production. Repair parts here refer to parts necessary for maintaining performance of the product.

# 10. Specifications

## Specification of the main unit

Model		GAS510C	
Online model		DL411C(recommended use), ADL312SC, GB211C	
Operating temperature		5-35°C (for indoor use only)	
Operating altitude		below an altitude of 2000m	
Operating humidity		≤75%RH	
Function		Provide closed anaerobic environment, condensate recovery solvent	
Safety device		Oxygen concentration meter, flammable gas alarm, spray nozzle detection switch, overheat prevention device	
Structure	System	Closed circulation type	
	Circulating gas	N <sub>2</sub>	
	Solvent recovery system	Condenser (evaporator) + refrigerator	
	Fine dust collecting system	Cartridge filter	
	Circulating blower	Roots blower	
	Refrigerator	Air cooling compressor	
		R404A 600g±10g	
	Solvent collecting container	2L flask	
	Meters	Filter differential pressure meter, flow meters (for introducing N <sub>2</sub> , for the oxygen sensor)	
	Filter	For protecting the O <sub>2</sub> sensor	
	Pump	For oxygen concentration measuring and circulation	
Compressor	For pushing the spray nozzle		
Control part	Evaporation temp. adjusting range	-20.0°C to 30.0°C	
	Evaporation temp. display	0.1°C	
	Circulating amount ※1	0.12~0.65 m <sup>3</sup> /min	
Spec.	External dimensions (mm) (WxDxH) ※2	710 x 950 x 1450	
	Power supply	200-230V~ 50/60Hz 5.5-12A	
	Weight	Approx. 240kg	
Accessories	<ul style="list-style-type: none"> <li>• Corrugated pipe 1300mm 1</li> <li>• Corrugated pipe 1000mm 1</li> <li>• Clamp (KF40) 4</li> <li>• Exhaust hose (gray) 1</li> <li>• Hoop (large) 1</li> <li>• Hose 5m 1</li> <li>• Hose 1.5m 2</li> <li>• Hoop (small) 6</li> <li>• Network cable (on-line) 1</li> </ul>	<ul style="list-style-type: none"> <li>• 2L collecting flask 1</li> <li>• Flask clamp 1</li> <li>• Connecting pipe (on-line) 1</li> <li>• Flat head screws (install the connecting pipe) 3</li> <li>• Warranty card 1</li> <li>• Instruction manual 1</li> <li>• Filter seal sticker 2</li> <li>• Sample hose Silicone: I.D. 2mmxO.D. 4mmx1m 2</li> <li style="padding-left: 20px;">Viton: I.D. 2mmxO.D. 4mmx1m 2</li> </ul>	


※1 Performance value: under the circumstance of AC230V power supply, room temperature 23°C ±5°C, humidity 65% RH± 5%, no load. The operating temperature range of this product is 5°C ~ 35°C.

※2 External dimensions do not include protruding parts.

# 11. Wiring diagram

## GAS510C


一次側 1/N/PE  
AC200~230V 50/60Hz



Symbol	Part Name	Symbol	Part Name
QR	Earth leakage breaker (30mA)	TB	Switching power
XD	Terminal block (switchboard)	BQ	Flammable gas alarm
2XD	Terminal block (back of unit body)	2BQ	Oxygen concentration meter
SF	Panel power switch	BT	Temperature sensor
2SF	Debug switch	PLC	Programmable logic controller
KFX	Solid state relay (refrigeration control)	HMI	Touch screen
2KFX	Solid state relay (remote start)	MAD	Analog input/output module (integrated)
KF	DC relay (air compressor)	AD	Analog input module (when used alone)
2KF	DC relay (Refrigeration unit)	DA	Analog output module (when used alone)
QCS	Roots blower controller (converter)	X	Network connector for spray dryer on-line
BE	Compressor current transmitter	KA	Combustible gas alarm
2BE	Temperature transmitter	2KA~7KA	Intermediate relay
GQ	air compressor	KH	Refrigeration control solenoid valve
2GQ	Refrigeration unit (air cooling)	2KH	AIR-IN solenoid valve
3GQ	Roots blower (blower)	3KH	Large flow N2 solenoid valve
4GQ	Vacuum pump for oxygen concentration detection	4KH	N2 normal open solenoid valve (small flow)
5GQ	Cooling fan	5KH	Spray gas source switching solenoid valve
6GQ	Cooling fan	6KH	Exhaust solenoid valve


# 12. System diagram

Piping system diagram



# 12. System diagram

Piping system diagram



# 12. System diagram

## List of materials

No.	Part name	Material	No.	Part name	Material
①	Condenser	SUS304	A	O-ring	Viton Adhesive: Silicone
②	Sensor	Stainless steel			
③	Ball valve	Stainless steel			
④	Clamp	POM	B	Packing	Silicone
⑤	Collecting flask	Glass	C	Hose	Teflon
⑥	Filter	Filter element	D	Spray nozzle	Label: Hycar (NBR) Teflon
		Filter case			
⑦	Differential pressure meter	ABS resin Diaphragm: Silicone Connecting hose: Viton	F	Aluminum honeycomb	Adhesive: phenol
⑧	Needle valve	Label assembly: Teflon Viton	G	Cap	Silicone
⑨	Flow meter ( for introducing N <sub>2</sub> )	Air contact assembly: Viton, Teflon Connecting hose : Teflon	H	Flowmeter	
⑩	Compressor	Cover case: PP Connecting hose : Teflon	I	Needle valve	Label assembly: Teflon Viton
⑪	Solenoid valve (for N <sub>2</sub> control)	Label assembly: Teflon	J	3-way valve ( for switching connection of GAS510C)	Label assembly: Teflon
⑫	Flow meter ( for measuring O <sub>2</sub> concentration)	Air contact assembly: Viton, Teflon Connecting hose: Viton	K	Solenoid valve (for pulse jet)	Label assembly: Teflon
⑬	Filter	Connecting hose: Viton	L	Packing	Silicone
⑭	Pump	Valve: Chloroprene Connecting hose: Viton			
⑮	O <sub>2</sub> sensor	Case: POM Packing: Viton Connecting hose: Viton			
⑯	Solenoid valve (for exhaust)	Label assembly: Teflon Connecting hose : Teflon			
⑰	Blower	Oil seal: NBR			
⑱	Solenoid valve ( for introducing N <sub>2</sub> )	Label assembly : Fluoro-rubber			
⑲	Solenoid valve (for suction)	Label assembly: Fluoro rubber			

## 13. Action principle

### Action principle

Please refer to P.53-55 「Piping system diagram」.

When the spray dryer is connected with the organic solvent recovery unit, a closed circulation system can be formed, and the organic solution can be safely treated after injection of inert gas.

The sample is fed from the appropriate container to the spray nozzle D by a liquid sending pump. In addition, the compressed air provided by the air compressor 10 and other equipment is regulated by the flow meter H and sent to the spray nozzle. At the tip of the nozzle, the compressed air is mixed with the sample, and the mixed sample is sprayed inside the drying chamber. At this time, the sample becomes droplet with a particle size of about 20 $\mu$ , and its surface area is 3,000cm<sup>2</sup> per 1L of sample.

On the other side, the air flow inside the unit can be recycled through the Roots fan in the organic solvent recovery unit. The air flow is heated to the set temperature through the heater to form a hot air, because the contact area of hot air and sample is very large, more than about 90% of the solvent is evaporated in the drying chamber instantly.

The sample that becomes fine particles by drying is sent to the cyclone part, where it is separated from the solvent vapor and collected into the product container A. The time from nozzle spraying to collecting into container is less than 0.5secs. Moreover, the sample particles are often surrounded by evaporated solvent vapor, and the temperature around the fine particles does not rise too high due to the vaporization heat. Therefore, substances that are not heat-resisting, such as enzyme, may also be powderized.

The evaporated solvent vapor, etc., through the circulating air stream, enters into the high-performance filter 6, which can filter the unrecycled sample particles. After the clean solvent vapor enters the condenser unit 1, the solvent vapor will liquefy rapidly at low temperature. Finally, collect the sample solvent into the waste liquid collecting flask 5.

After the sample is treated, there is also solvent residue in the circulating pipeline system. Please run for a longer period of time to ensure solvent recovery. Then, start the air suction button, turn on the air suction solenoid valve 19, when the internal oxygen concentration returns to more than 18%, it's able to stop the unit.

The temperature conditions in the experiment will be set and displayed on the display panel of spray dryer through the inlet temperature sensor and the outlet temperature sensor.

The condensation temperature conditions in the experiment will be set and displayed on the display panel of organic solvent recovery unit through the condensation temperature sensor.

The oxygen concentration conditions in the experiment will be displayed on the display panel of organic solvent recovery unit through the oxygen concentration analyzer sensor.

In addition, you can also select the optional circulating air flowmeter, so that the air flow of drying sample will also be displayed on the display panel of organic solvent recovery unit.

## 14. Replacement parts table

	Part name	Specification	Manufacturer	Code No
※	Filter	GAS41-40040 for refrigerator	YSC	B040300003
	Flammable gas detector	TJ-PLT119-EX	YSJ	A020299002
	Axial flow fan	SJ1238HA2BAL for thermal discharge	YSJ	A080104012
	Temperature sensor	GAS410C_03_01-02	YSJ	H100101028
	Roots blower	IRS-32A	YSJ	B040201001
	Miniature vacuum pump	D35S-41J-0000 for oxygen sensor	YSJ	A041700055
※	Activated carbon filter set	GAS41-40570 for oxygen sensor	YSC	B081601001
	Flow meter for introducing N <sub>2</sub>	LZB-07A10MT 5-50L/min	YSJ	A040409005
	Flow meter for oxygen sensor	LZM-6T 0.1-1.5L/min	YSJ	A040499023
	Solenoid valve	CKD AB41-04-8-M-AC100V for air suction	YSJ	A040403015
	Solenoid valve	CKD AB41-03-7-F-AC100V for air exhaust	YSJ	A040403014
	Solenoid valve	CKD AB31-01-2-M-AC100V for introducing N <sub>2</sub>	YSJ	A040403013
	Solenoid valve	CKD AB41-02-5-F-AC100V for controlling N <sub>2</sub>	YSJ	A040403004
	Solenoid valve	1328N/2S030 For refrigeration control	YSJ	A031800003
	Differential pressure meter	DG87-641-1C 0~1kPa	YSC	B042300001
※	Filter element	RE-205-90-FB 0.3μ	YSC	B040300001
※	Fluoro rubber hose	φ4×φ6 specify length	YSC	B080807049
※	Teflon hose	φ6×φ8 specify length	YSJ	A080807007
※	Teflon hose	φ10×φ12 specify length	YSJ	A080807006
	Corrugated pipe	40KF L=160	YSJ	A041500019
	Corrugated pipe	40KF L=1500	YSJ	A041500006
	Corrugated pipe	40KF L=1000	YSJ	A041500003
	Corrugated pipe	40KF L=500	YSJ	A041500008
	Corrugated pipe	40KF L=1300 for external connection	YSJ	A041500004
	Corrugated pipe	40KF L=1000 for external connection	YSJ	A041500003
	Clamp	40KF	YSJ	A041500078
※	Center ring	40KF	YSJ	A041500077
	Duct hose	φ25×2m	YSC	B080807030
※	Hose	Φ6×Φ11×5m	YSC	B080807050
※	Liquid sending hose	φ 2×φ4(PTFE)	YSJ	B080913005
※	Liquid sending hose	φ 2×φ4×1m	YSJ	B080807050

## 14. Replacement parts table

Part name	Specification	Manufacturer	Code №
PLC	CP2E-N30DR-A	YSJ	A020300069
Analog input/output unit (integrated)	CP1W-MAD44	YSJ	A020399053
Analog input unit (alone)	CP1W-AD042	YSJ	A020399047
Analog output unit (alone)	CP1W-DA042	YSJ	A020399065
Touch screen	NB7W-TW11B	YSJ	A020400014
CAT6 gigabit network cable	2m	YSJ	A120102034
CAT6 gigabit network cable	1m	YSJ	A120102032
CAT6 gigabit network cable	0.5m	YSJ	A120102031
Switching power	LRS-100-24	YSJ	A010801045
Intelligent temperature transmitter	NHR-213	YSJ	A010599023
Frequency converter	3G3JZ-AB007	YSJ	A020501001
Stationary oxygen detector	HT-FX100	YSJ	A020799033
DC relay	HF13F/024Z21D	YSJ	A011001011
DC relay	HF41F/12-ZS+41F-1Z-C2-1	YSJ	A011001027
DC relay	HF41F/24-ZS+41F-1Z-C2-1	YSJ	A011001028
Solid state relay	SSK4A2032+ guide rail base	YSJ	A011006029
※ Filter seal sticker	130*130*0.5 (thickness)	YSJ	A089900040
Earth leakage circuit breaker	BV-DN IP+N 16A 30mA	YSJ	A010410004
Air cooling condensing unit	CAJZ2446ZBR	YSJ	A030101013
Oil-free air compressor	688CGHI44	YSC	A040201008
Axial flow fan	SJ2206LA2BAT for heat dissipation	YSJ	A080104030

Note: Parts marked with ※ are consumable parts.

# 15. List of Dangerous Substances



Never use explosive substances, flammable substances and substances that include explosive or flammable ingredients in this unit. Otherwise explosion or fire may result

Explosive substance	Explosive substance	① Nitroglycol, glycerine trinitrate, cellulose nitrate and other explosive nitrate esters
		② Trinitrobenzen, trinitrotoluene, picric acid and other explosive nitro compounds
		③ Acetyl hydroperoxide, methyl ethyl ketone peroxide, benzoyl peroxide and other organic peroxides
Flammable substances	Explosive substances	Metal "lithium", metal "potassium", metal "natrium", yellow phosphorus, phosphorus sulfide, red phosphorus, celluloids, calcium carbide (a.k.a, carbide), lime phosphide, magnesium powder, aluminum powder, metal powder other than magnesium and aluminum powder, sodium dithionous acid (a.k.a., hydrosulphite)
	Oxidizing substances	① Potassium chlorate, sodium chlorate, ammonium chlorate, and other chlorates
		② Potassium perchlorate, sodium perchlorate, ammonium perchlorate, and other perchlorates
		③ Potassium peroxide, sodium peroxide, barium peroxide, and other inorganic peroxides
		④ Potassium nitrate, sodium nitrate, ammonium nitrate, and other nitrates
		⑤ Sodium chlorite and other chlorites
		⑥ Calcium hypochlorite and other hypochlorites
	Flammable substances	① Ethyl ether, gasoline, acetaldehyde, propylene chloride, carbon disulfide, and other substances with ignition point at a degree 30 or more degrees below zero.
		② n-hexane, ethylene oxide, acetone, benzene, methyl ethyl ketone and other substances with ignition point between 30 degrees below zero and less than zero.
		③ Methanol, ethanol, xylene, pentyl acetate, (a.k.a.amyl acetate) and other substances with ignition point between zero and less than 30 degrees.
④ Kerosene, light oil, terebinth oil, isopentyl alcohol(a.k.a. isoamyl alcohol), acetic acid and other substances with ignition point between 30 degrees and less than 65 degrees.		
Combustible gas	Hydrogen, acetylene, ethylene, methane, ethane, propane, butane and other gases combustible at 15°C at one air pressure.	

(Quoted from the separate table 1 in Article 6, the enforcement order of the Industrial Safety and Health Law)

## Responsibility

Please follow the instructions in this document when using this unit. Yamato Scientific has no responsibility for the accidents or breakdown of device if it is used with a failure to comply. Never conduct what this document forbids. Unexpected accidents or breakdown may result in.

## Note

- ◆ The contents of this document may be changed in future without notice.
- ◆ Any books with missing pages or disorderly binding may be replaced.

Instruction Manual  
Solvent Recovery Unit  
GAS510C  
Version 1      Dec. 06, 2024  
Revision      Apr. 07, 2025

---

YAMATO SCIENTIFIC CO., LTD.  
Harumi Triton Square Y-36F, 1-8-11 Harumi,  
Chuo-ku, Tokyo 104-6136, Japan  
Tel : +81-3-5548-7122  
Fax : +81-3-5548-0132  
<https://www.yamato-scientific.com/>